


NEW

PopLok™

2.8 & 3.3mm All PEEK Knotless Anchor



- Can be Single loaded with #2 Hi-Fi® Suture or Double loaded with #0 Hi-Fi Suture
- Combination press-fit rib design and subcortical wings enhance fixation
- Sutures can be individually tensioned separately from anchor seating
- Anchor is cannulated to channel growth factors

Clinical Use Hip Labral Repair

CONMED
LIMVATEC

Unique hole-finding tip

Easy-to-load Suture Tab

Ideal length for hip repair

Ergonomic PopLok Handle Design

2.8 & 3.3mm PopLok™ Suture Anchors:

These all-PEEK knotless anchors are constructed utilizing the popular PopLok design and are ideal for Labral Repair.

The 2.8mm design boasts the smallest size double loaded knotless anchor available on the market!



Actual Size: 2.8mm

Mean Pull-Out Strengths

2.8mm – 174 Newtons

3.3mm – 247 Newtons

ORDERING INFORMATION

PopLok™ Suture Anchor

- GKP-2801** 2.8mm PopLok suture anchor with one #2 Hi-Fi® suture
- GKP-2802** 2.8mm PopLok suture anchor with two #0 Hi-Fi sutures
- GKP-3301** 3.3mm PopLok suture anchor with one #2 Hi-Fi suture
- GKP-3302** 3.3mm PopLok suture anchor with two #0 Hi-Fi sutures

PopLok Instrumentation

- BKL-00M** 2.8/3.3mm PopLok Drill Bit
- BGU-00M** PopLok Guide



CrossFit™ and PopLok™ are comprised of PEEK-OPTIMA® polymer from Imbico® Biomaterial Solutions. PEEK-OPTIMA and IMBICO are registered trademarks of Imbico Ltd. All Rights are Reserved.



This product contains Dynema® Purity. Dynema® Purity is a registered trademark of Royal DSM N.V.

Visit Linvatec.com today to learn more.

CONMED CORPORATION PRODUCT AREAS:

ARTHRISCOPY • ELECTROSURGERY • ENDOSCOPY • ENDOSURGERY • GASTROENTEROLOGY • INTEGRATED SYSTEMS • PATIENT CARE • POWERED INSTRUMENTS • PULMONOLOGY



11311 Concept Boulevard
 Largo, FL 33773-4908
 (727) 392-6464
 Customer Service: 1-800-237-0169
 FAX: (727) 399-5256
 International FAX: +1 (727) 397-4540
 email: customer_service@linvatec.com
www.linvatec.com